



Step by Step Guide

SET UP


Install and Set up




Loading screen




Select Child/Person with ASD avatar.




Sing in / Sing Up methods.



Create a new passcode to protect the menu in Child mode.




Select your preferred language.






Select Parent/caregiver avatar.




Child role level selection.


Caregiver mode. Speech to image.



Parent mode. Tap the "Record" button and record your sentence.




Tap the "Stop" button and wait - images will appear.




The sentence is simplified by removing words that cannot be represented with images.

Tap "Clear" to restart the process.


Caregiver mode. Missing image - Missing image generator.




Parent mode. Tap the "Record" button and record your sentence.




Select the most relevant image.




Tap the "Stop" button and wait - images will appear.




Choose the most relevant category and subcategory.




The sentence is simplified by removing words that cannot be represented with images. One image is missing in our library. Tap the "Missing image" icon.



After you save the image with the assigned categories, it will become public.




Tap the "Generate" button to create four images with AI assistance.




Instead of the "Missing image" icon, you will now see the image.


Caregiver mode. Missing image - Own image upload.




Parent mode. Tap the "Record" button and record your sentence.




Tap the "Stop" button and wait - images will appear.




The sentence is simplified by removing words that cannot be represented with images. One image is missing in our library. Tap the "Missing image" icon.




Tap the "Generate" button to create four images with AI assistance.




If you do not like images, you can upload your own image for this word.




Choose whether you want to upload an image from your files or capture an image with your camera.



Choose the most relevant category and subcategory.




After you save the image with the assigned categories, it will become private.




Instead of the "Missing image" icon, you will now see the image.


Caregiver mode. Missing image flow - Draw pad.




Parent mode. Tap the "Record" button and record your sentence.




Draw the image.




Tap the "Stop" button and wait - images will appear.




Tap the "Use" button.




The sentence is simplified by removing words that cannot be represented with images. One image is missing in our library. Tap the "Missing image" icon.




Choose the most relevant category and subcategory.




Tap the "Generate" button to create four images with AI assistance.



After you save the image with the assigned categories, it will become private.




If you don't like the images, you can use the drawing pad to create your own image for this word.




Instead of the "Missing image" icon, you will now see the image.


Caregiver mode. Missing image flow - Draw pad - Drawing improvement with Al.




Parent mode. Tap the "Record" button and record your sentence.




Tap the "Stop" button and wait - images will appear.




The sentence is simplified by removing words that cannot be represented with images. One image is missing in our library. Tap the "Missing image" icon.




Tap the "Generate" button to create four images with Al assistance.




If you don't like the images, you can use the drawing pad to create your own image for this word.




Draw the image.




Tap the "Upgrade" button.




Tap the "Generate" button to create four images based on your drawing with AI assistance.



Choose the most relevant image.



Choose the most relevant category and subcategory.



After you save the image with the assigned categories, it will become private.




Instead of the "Missing image" icon, you will now see the image.

CHILD MODE

Individual with ASD mode. Image to speech / Image to visual sentence.

Beginner level. Only Favourite images available - Favourite images selected by parent.



Child mode: Beginner level. Empty image placeholders will appear until the parent adds images for use (refer to the "Favourite Images" instructions).



Beginner level. Child can use images to provide answer or compose sentences.




The child must tap on the relevant image to compose answers.


CHILD MODE

Individual with ASD mode. Image to speech / Image to visual sentence.

Advanced level. For experienced user. Unlimited image library available. Favourite images selected by parent and placed to Favourites category.




Child mode: Advanced level. The child can use the image library, and categories can help find images. The 'Favourite' category allows the parent to select images for quick use (refer to the 'Favourite Images' instructions). Tap to choose a category.



Some categories have sub-categories. Tap to choose a sub-category.




Go Back to choose a different category and sub-category.




Compose a sentence using images. Tap "Clear" to remove the sentence. Tap on an image (at the top of the screen) to remove it.

PRIVATE IMAGE UPLOAD

Parent mode. Replace a public image with a private image.




Parent mode. Tap the "Record" button and record your sentence. Tap the "Stop" button and wait - images will appear. The sentence is simplified by removing words that cannot be represented with images.




If you want to replace a public image (from our public image library) with your own (private image), tap on the image.


Before




Choose whether you want to upload an image from your files or capture an image with your camera.



Tap "Save" and the public image will be replaced, but it will remain private (only you will be able to see it).




You can find your uploaded images later in the Menu. "Uploaded images" section.




After

PRIVATE IMAGE - REMOVE OR REPLACE


Menu section: Uploaded images. Replace or remove private images.



Go to the menu and select "Uploaded Images".




You can view all private images you've uploaded using the camera, your smartphone files, or the drawing pad.




You can remove your private image and you will see again public image for this word.

Before




Choose whether you want to upload an image from your files or capture an image with your camera.



Tap "Save" and your private image will be replaced, and it will remain private (only you will be able to see it).




Now, in "Uploaded Images", you will see the updated image.




After

FAVOURITE IMAGES


Menu section: Manage Favourites. Favourite images selected by parent.




Go to the menu and select "Manage Favourites".




To add new images, tap "Add New".



Search for the word you want to add to the Favourite folder.




Press "Add" if you are happy with the result, or press "Clear" to start again.




After you add the image, it will appear in "Manage Favourites".

Beginner level



Beginner level: Favourited images will appear on the main screen.

Advanced level




Advanced level: Favourited images will appear in the "Favourites" category. Tap on it to view the images.




"Favourites" category.

SWITCH MODE

Menu section: Switch mode. Switch role. Select level. Quick role switch toggle ON or OFF.



Go to the menu and select "Settings".



Role Parent Child Quick role switch toggle Toggle

two devices.

← Switch mode

Same level will also be applied to the

9

Level


Beginner

Advanced


connected device.

You can choose the Level and Mode. To quickly switch roles from the home screen (parent to child; child to parent), turn the toggle On.


It makes easy to switch between two roles, if you use one device. Uncheck, if you use



Toggle is On. Use the Toggle to switch the role.



Toggle is On. Use the Toggle to switch the role.




 \equiv


Toggle is Off.

SETTINGS

Menu section: Settings. Notifications. Voice accompaniment. Language. Avatars. Child lock passcode.



Go to the menu and select "Settings".



- Turn app notifications On or Off.
- Choose whether you want to hear word pronunciation when tapping on an image in Child mode (Image to Speech function).
- Select your preferred language.
- Select avatars for Child or Parent.
- Create a new passcode to protect the menu in Child mode (this works only if the Toggle is Off in "Quick Role Switch" under Menu > Switch Mode).